

SINGLE CHANNEL PIPETTES

A completes our range of micropipettes with the Fully Autoclavable feature while maintaining the light weight and ergonomic design and ensuring high accuracy and precision.

FEATURES



Fully Autoclavable

Click-stop digital system for volume setting.

Ergonomic Design ensuring light weight & soft plunger movement.

High accuracy and precision guaranteed

Two step plunger operation allows "reverse pipetting technique".

Individual serial No. permanently marked on each pipette.

Calibration as per DIN 12650 & ISO 8655 standards.

Calibration report enclosed with every pipette.

Color coded for easy identification.

"CE" Mark

MULTICHANNEL PIPETTES

- Soft force / leak free tip loading and ejection.
- Easily removable & fully autoclavable manifold.
- Easy Right- or left-hand operation because of 360° rotatable Manifold.
- One hand tip ejection.
- Attractive color coding



VARIABLE VOLUME PIPETTE

Cat. No.	Volume Range [μl]	Increments [μl]	Test Volume [μl]	Inaccuracy [μl] %	Imprecision [μl] %	Price
BS-125	0.1 - 2.5 μl	0.0020	0.25	12	6	P.O. R
			2.5	2.5	1.6	
BS-100	0.5-10 μl	0.0200	1	2.5	1.8	P.O. R
			10	1	0.6	
BS-200	5-50 μl	0.1000	5	2	2	P.O. R
			50	0.6	0.3	
BS-500	10-100 μl	0.1000	10	3	1	P.O. R
			100	0.8	0.2	
BS-600	100-1000 μl	0.1000	100	3	0.6	P.O. R
			1000	0.6	0.2	
BS-700	2-20 μl	0.0200	2	3	2	P.O. R
			20	0.9	0.4	
BS-800	20-200 μl	0.2000	20	2.5	0.7	P.O. R
			200	0.6	0.2	
BS-905	500-5000 μl	10.0000	500	2.4	0.6	P.O. R
			2500	1.2	0.25	
			5000	0.6	0.15	
BS-1000	1000-10000 μl	20.0000	1000	3	0.6	P.O. R
			5000	1.2	0.3	
			10000	0.6	0.2	

FIXED VOLUME PIPETTE

Cat. No.	Volume Range [μl]	Test Volume [μl]	Inaccuracy [μl] %	Imprecision [μl] %	Price
BS-2.5	2.5 μl	2.5	2	1.6	P.O. R
BS-0005	5 μl	5	1.3	1.2	P.O. R
BS-0010	10 μl	10	1.2	0.8	P.O. R
BS-0020	20 μl	20	1	0.5	P.O. R
BS-0025	25 μl	25	1	0.3	P.O. R
BS-0050	50 μl	50	0.7	0.3	P.O. R
BS-0100	100 μl	100	0.6	0.2	P.O. R
BS-0200	200 μl	200	0.6	0.2	P.O. R
BS-0250	250 μl	250	0.6	0.2	P.O. R
BS-0500	500 μl	500	0.6	0.2	P.O. R
BS-1000	1000 μl	1000	0.6	0.2	P.O. R
BS-2000	2000 μl	2000	0.3	0.15	P.O. R
BS-5000	5000 μl	5000	0.3	0.15	P.O. R
BS-10000	10000 μl	10000	0.6	0.2	P.O. R

MULTICHANNEL PIPETTE

Cat. No.	8 & 12 Channel Volume Range [μl]	Increments [μl]	Test Volume [μl]	Inaccuracy [μl] %	Imprecision [μl] %	Price
BS8-010 BS12-010	0.5-10 μl	0.0200	10	2	1	P.O. R
			5	4	2	
			1	8	5	
BS8-020 BS12-020	2 -20 μl	0.0200	20	2	1.6	P.O. R
			10	3	2	
			2	7	3	
BS8-050 BS12-050	5-50 μl	0.1000	50	1	0.7	P.O. R
			25	1.7	1	
			5	3	2	
BS8-100 BS12-100	10-100 μl	0.1000	100	0.8	0.3	P.O. R
			50	1	0.8	
			10	3	2	
BS8-200 BS12-200	20-200 μl	0.2000	200	0.7	0.25	P.O. R
			100	1	0.4	
			20	5	1.4	
BS8-300 BS12-300	30-300 μl	0.2000	300	0.6	0.3	P.O. R
			150	1	0.5	
			30	3	1	